



ANESTHESIA

RETROBULBAR, Atkinson Tip





Code	Description	Contents
111638	23G x 1½ in. (0.6 x 38mm)	10/Box
111637	25G x 1½ in. (0.5 x 38mm)	10/Box

RETROBULBAR, Sharp Tip





Code	Description	Contents
111278	25G x 1 in. (0.5 x 25mm)	10/Box
111282	30G x ½ in. (0.3 x 12.5mm)	10/Box

PERIBULBAR, Atkinson Tip



Code	Description	Contents
115144	23G x 7/8 in. (0.6 x 22mm)	10/Box
115027	25G x 7/8 in. (0.5 x 22mm)	10/Box
111642	25G x 1¼ in. (0.5 x 32mm)	10/Box



SUB-TENON



Code	Description	Contents
125176	19G x 7/8 in. (1.1 x 22mm), open tip	10/Box



ANTERIOR CHAMBER IRRIGATION

ANGLED TIP, 4mm (Rycroft)





Code	Description	Contents
121281	20G (0.9mm)	10/Box
121639	25G (0.5mm)	10/Box
101280	27G (0.4mm)	10/Box
121273	30G (0.3mm)	10/Box

ANGLED TIP, 6mm (Rycroft)





Code	Description	Contents
121287	20G (0.9mm)	10/Box
121286	27G (0.4mm)	10/Box
121276	30G (0.3mm)	10/Box

ANGLED TIP, 7mm (Rycroft)




Code	Description	Contents
121637	25G (0.5mm)	10/Box
121290	27G (0.4mm)	10/Box
121277	30G (0.3mm)	10/Box

ANGLED TIP, 8mm (LASIK)



Code	Description	Contents
121293	27G (0.4mm)	10/Box



ANGLED TIP, 10mm (LASIK)



Code	Description	Contents
121280	25G (0.5mm)	10/Box
121274	27G (0.4mm)	10/Box


CAPSULE POLISHING

SILICONE TIP, Angled 7mm




Code	Description	Contents
175041	23G (0.6mm)	10/Box
175005	27G (0.4mm)	10/Box

POSTERIOR, Angled 9mm, Top Port





Code	Description	Contents
171604	21G (0.8mm)	10/Box

POSTERIOR, Angled 7mm, Open Tip



Code	Description	Contents
171614	21G (0.8mm)	10/Box


OLIVE TIP



Code	Description	Contents
171589	23G (0.6mm)	10/Box

CORTEX REMOVAL


CORTEX ASPIRATOR (Binkhorst)



Code	Description	Contents
121295	23G (0.6mm), left	10/Box
121294	23G (0.6mm), right	10/Box


CAPSULOTOMY

FORMED





Code	Description	Contents
131599	23G x 5/8 in. (0.6 x 16mm)	10/Box
131610	25G x 5/8 in. (0.5 x 16mm)	10/Box
131618	27G x 5/8 in. (0.4 x 16mm)	10/Box

ANGLED 45°




Code	Description	Contents
135049	25G x 5/8 in. (0.5 x 16mm)	10/Box
135093	27G x 5/8 in. (0.4 x 16mm)	10/Box

REVERSE BEND




Code	Description	Contents
135115	25G x 5/8 in. (0.5 x 16mm)	10/Box

SHORT





Code	Description	Contents
135004	25G x 5/8 in. (0.5 x 16mm)	10/Box
135008	27G x 5/8 in. (0.4 x 16mm)	10/Box

SIDE-CUT



Code	Description	Contents
131602	25G x 5/8 in. (0.5 x 16mm)	10/Box
131609	27G x 5/8 in. (0.4 x 16mm)	10/Box



PUSH-PULL, Formed



Code	Description	Contents
135042	27G x 5/8 in. (0.4 x 16mm)	10/Box



HYDRODISSECTION

FLAT TIP, Angled 7mm





Code	Description	Contents
125024	23G (0.6mm)	10/Box
125036	25G (0.5mm)	10/Box
125162	27G (0.4mm)	10/Box

FLAT TIP, Angled 8mm





Code	Description	Contents
125037	25G (0.5mm)	10/Box
125158	27G (0.4mm)	10/Box
121592	30G (0.3mm)	10/Box

HORIZONTAL FLAT TIP, Curved




Code	Description	Contents
125155	27G (0.4mm)	10/Box

VERTICAL FLAT TIP, Curved



Code	Description	Contents
125099	27G (0.4mm)	10/Box



J-SHAPED, Vertical Flat Tip



Code	Description	Contents
125029	27G (0.4mm)	10/Box



HYDRODISSECTION

ANGLED STANDARD TIP (Chang)





Code	Description	Contents
121198	27G (0.4mm), 1.5mm tip	10/Box

CORTICAL CLEAVING, Flat Tip, Angled 11mm





Code	Description	Contents
105158	27G (0.4mm)	10/Box

CORTICAL CLEAVING, Round Tip, Angled 11mm





Code	Description	Contents
125157	27G (0.4mm)	10/Box

CORTICAL CLEAVING, Round Tip, Angled 12mm



Code	Description	Contents
125197	30G (0.3mm)	10/Box



TAPERED, Angled 1mm



Code	Description	Contents
125159	25 to 30G (0.5 to 0.3mm)	10/Box



IRRIGATION AND AIR INJECTION

VISCOELASTIC, Angled





Code	Description	Contents
125009	23G x $\frac{7}{8}$ in. (0.6 x 22mm)	10/Box
125026	27G x $\frac{7}{8}$ in. (0.4 x 22mm)	10/Box

VISCOELASTIC, Formed





Code	Description	Contents
125007	27G x $\frac{7}{8}$ in. (0.4 x 22mm)	10/Box

WASHOUT, Angled 8mm (Bishop-Harmon)



Code	Description	Contents
165040	19G x $\frac{5}{16}$ in. (1.1 x 8mm)	10/Box

OLIVE TIP, Curved



Code	Description	Contents
165147	23G x $\frac{7}{8}$ in. (0.6 x 22mm)	10/Box

NUCLEUS REMOVAL



VECTIS, 1 Port



Code	Description	Contents
141619	24G (0.55mm)	10/Box

IRRIGATION AND ASPIRATION



SIMCOE, Standard



Code	Description	Contents
161571	23G (0.6mm), disposable	10/Box
161573	23G (0.6mm), reusable	1/Box

ANTERIOR CHAMBER MAINTENANCE


RIBBED OPEN TIP



Code	Description	Contents
125061	20G (0.9mm)	3/Box
125062	23G (0.6mm)	3/Box



OCULOPLASTIC

LACRIMAL, Curved



Code	Description	Contents
181615	21G x 1 ¹ / ₈ in. (0.8 x 28mm)	10/Box
181276	26G x 1 ¹ / ₈ in. (0.45 x 28mm)	10/Box



LACRIMAL, Shahinian



Code	Description	Contents
185033	23G (0.6mm), bullet tip with top port	5/Box



REFRACTIVE

BURATTO, Formed 9mm




Code	Description	Contents
205242	25G (0.5mm), 2 ports, open round tip	10/Box

GUELL, Formed 9mm




Code	Description	Contents
205240	25G (0.5mm), 6 ports, open round tip	10/Box

HERSH, Formed 9mm



Code	Description	Contents
205252	23G (0.6mm), open spatula tip	10/Box
205251	25G (0.5mm), open spatula tip	10/Box


LANE, Angled



Code	Description	Contents
205248	25G (0.5mm), 2 ports, open spatula tip	10/Box



REFRACTIVE

LINDSTROM, Formed 9mm





Code	Description	Contents
205250	25G (0.5mm), 6 ports, closed tip	10/Box

NICHAMIN, Formed 9mm




Code	Description	Contents
205254	25G (0.5mm), horizontal open flat tip	10/Box

QUAD PORTED, Curved





Code	Description	Contents
205249	23G (0.6mm), 4 ports, closed tip	10/Box

SLADE, Angled





Code	Description	Contents
205246	25G (0.5mm), open spatula tip	10/Box

SOLOMON, Formed 9mm



Code	Description	Contents
205280	25G (0.5mm), 6 ports, open spatula tip	10/Box


VIDAURRI, Angled



Code	Description	Contents
205256	25G (0.5mm), 2 angled tip extensions, 8 ports, closed tips	10/Box


VITREORETINAL

PERFLUOROCARBON




Code	Description	Contents
184923	30G (0.3mm) inner/23G (0.6mm) outer	10/Box

SOFT TIP, Straight



Code	Description	Contents
185185	20G (0.9mm), 2mm extension	10/Box
185195	23G (0.6mm), 2mm extension	10/Box

VITREORETINAL BRUSH



Code	Description	Contents
185138	20G (0.9mm), 1mm extension	10/Box
185141	23G (0.6mm), 1mm extension	10/Box

LABELING

STERILE R

Sterilization Method: "R" indicates sterilization by gamma irradiation; "EO" indicates ethylene oxide sterilization



Single-use symbol



Indicates year/month of product expiration
Example: 2012-01 = January, 2012

LOT

Manufacturing lot for traceability



See instructions for use













European Union CE mark and Notified Body ID number



Recycle symbol

HUBS

Surgistar hubs are color-coded in accordance with ISO 6009 standards, with the exception of 30G (clear).

Color		Gauge	mm
Clear		30G	0.30
Grey		27G	0.40
Brown		26G	0.45
Orange		25G	0.50
Purple		24G	0.55
Blue		23G	0.60
Black		22G	0.70
Green		21G	0.80
Yellow		20G	0.90
Cream		19G	1.10