

# 2022 STERILISATION TRAYS CATALOGUE

**SurgiVision**  
SurgiVision Pty Ltd



*Call 1300 995 132 now to order  
OR email [sales@surgivision.com.au](mailto:sales@surgivision.com.au)*

[www.surgivision.com.au](http://www.surgivision.com.au)

# Micro Ophthalmic Baskets

## Micro Ophthalmic Baskets

Micro ophthalmic baskets are specially designed for delicate micro instruments.

ITEM CODE	SIZES (mm)
SV90-22000	240x250x38
SV90-22001	480x250x38
SV90-22002	400x250x38



- With silicon inserts on bottom of baskets to accommodate micro instruments.
- Lids are provided with silicon strips to grip the instruments.
- Side walls are made of perforated sheet.
- User friendly locking handle.

# Micro Instrument Cassettes

## Spring Holding Cassette

Specially designed to minimise the contact area while locking of instruments inside the cassette. Each instrument is held by separate spring. Further it can equally hold thicker and thinner instruments in one arrangement of the cassette.

ITEM CODE	SIZES (mm)
SV-4242	178x64x32
SV-4243	178x114x32
SV-4244	203x140x38
SV-4245	203x280x38
SV-4246	220x125x40
SV-4247	250x220x40



Call 1300 995 132 now to order  
OR email [sales@surgivision.com.au](mailto:sales@surgivision.com.au)

[www.surgivision.com.au](http://www.surgivision.com.au)



# Micro Instrument Cassettes


## Fine Mesh Cassette

Fine Mesh cassettes are ideal when instruments are required to be cleaned by all around water flow.


ITEM CODE	SIZES (mm)
SV-4148	178x64x32
SV-4149	178x114x32
SV-4150	203x140x38




ITEM CODE	SIZES (mm)
SV-3010	180x60x22
SV-3011	180x125x22



ITEM CODE	SIZES (mm)
SV-3012	184x145x34
SV-3013	184x74x34
SV-3014	184x74x34



ITEM CODE	SIZES (mm)
SV-3017	280x183x32



Call 1300 995 132 now to order  
OR email [sales@surgivision.com.au](mailto:sales@surgivision.com.au)

[www.surgivision.com.au](http://www.surgivision.com.au)


# Instrument Cassettes

■ **ITEM CODE**      **SIZES (mm)**  
 SV-3016      ▶      280x183x32

For 16 instruments with partition plate

■ **ITEM CODE**      **SIZES (mm)**  
 SV-3017      ▶      280x183x32

For 20 instruments



■ **ITEM CODE**      **SIZES (mm)**  
 SV-3018      ▶      280x205x32



■ **ITEM CODE**      **SIZES (mm)**  
 SV-3019      ▶      280x183x32



Silicon inserts available in different colours.



Green



Orange



Yellow



Red



Blue