




BALANCED SALT SOLUTION

Intraocular irrigation solution



FOR OPHTHALMIC PROCEDURES

Eye-Lotion

BALANCED SALT SOLUTION



Distributed by

SurgiVision

SurgiVision Pty Ltd

Call 1300 995 132 now to order
OR email sales@surgivision.com.au

www.surgivision.com.au

CE0483

Balanced Salt Solution

Sterile physiological solution for ophthalmic use.

Intended Purpose

To irrigate the eye during ophthalmic procedures that require irrigation or moisture retention.

The period of use is limited by the duration of the procedure.

Please observe the instructions for use.

1 ml solution contains

Sodium Chloride	6.40 mg
Calcium Chloride 2H ₂ O	0.48 mg
Sodium Acetate 3H ₂ O	3.90 mg
Potassium Chloride	0.75 mg
Magnesium Chloride 6H ₂ O	0.30 mg
Sodium Citrate 2H ₂ O	1.70 mg

Type	Pack Sizes	Product Code
Plastic-Bottle	10 x 250ml	SW 1394
Plastic-Bottle	10 x 500ml	SW 1395
Plastic-Bottle	50 x 20ml	SW 1396



500 ml

SW 1395



250 ml

SW 1394



20 ml

SW 1396

Certified to EN ISO 13485:2016

TGA ARTG N° 409903

Call 1300 995 132 now to order
OR email sales@surgivision.com.au
www.surgivision.com.au

SurgiVision
SurgiVision Pty Ltd

serumwerk
bernburg